

This Procedure is for those who wants to use eeglab to process the EGI data.

Tested software version: Matlab 2015b, Nestation 5.3.0.1, EEGLAB V13.6.5b


Before we export data to EEGLAB, we need to know some difference to treat Events name between Netstation and EEGLAB.

Topic 1 Markup function

Net Station events contain key lists, in other words, mini databases. Using a simple standard/target experiment as an example, in Net Station all stimuli events might be named “stim.” The distinction between standards and targets can’t be determined from the name of the event. It can only be determined from the key list in the events.

In most systems, including FieldTrip, events only have names. They don’t have key lists. To make the jump from key-list events to non-keylist events, you must use Net Station’s segmentation markup. Segmentation markup adds events to the recording that you can use to distinguish different conditions. To use segmentation markup, create a segmentation specification, and check the “Mark Up File” checkbox in the segmentation specification editor. When you run segmentation using this specification, instead of segmenting the file, new events will be added. For example, if all your standard and target stimuli are named “stim”, you would create a segmentation specification just as you would for segmenting this file into standard and target categories. Then, if you check the “Mark Up File” checkbox, the specification will cause new events to be added to the file instead of segmentation. Then, you can add events called “stnd” for standard, and “targ” for target.

Note: Use 4 ASCII characters to name the events. Do not include Space.




Topic 2: Load location file and Average reference in EEGLAB

Then you need to consider if you want to make use of the “bad channel replacement” tool in NetStation. Sometimes it’s very useful, especially when there are some bad channels while you do want to retain them.

If you do not want it and just need to finish all the work in EEGLAB, please go to Topic 2 . If you think the toos is necessary for you, then go to Topic 3.

Step 1: export your file to Simple binary file. Launch Netstation tools, and create a file export tool.


You can check the options what you want to keep. For example, if you want to see Ref channels in EEGLAB, then check option"Export Reference Channel". Then name the tool, and click "Create" button.



Step 2: Run the tool on the file you want to convert. After that, you will see a file whose format is raw.


Note:

- EEGLAB does not recognize breaks in recording (Epoch breaks) when reading .raw files.
- To avoid recording breaks with EPrime for experiment control, add "false" to the end of parameter line in NSInit object, such as following:
c, "on", CellList, "socket", "10.10.10.42", false
And do not put NSSStopRecording object except at the end of experiment.




Step 3: Launch eeglab and import the data.

Go to File/import data/Using EEGLAB functions and plugins/From nestation binary simple file:



The following windows will tell you EEGLAB will use the default location file, just click Ok. Then name your file:



Step 4: Load the location file.

Please visit www.egi.com and enter knowledge center by clicking "Training and Technical support" tab. Then search the keyword "location". You will see a list of topics related to keyword "location". Click the link named **"3D Sensor Position files for visualizing data in third party software"**

The screenshot shows the EGI Knowledge Center website. The top navigation bar includes links for Home, Contact, Education, and Training and Technical Support. Below this is a secondary navigation bar with links for Company, Clinical Division, Research Division, Neuromodulation, Science at EGI, and Investors. A search bar is located on the right side of this bar, containing the text 'location'. The main content area features the EGI logo and the text 'ELECTRICAL GEODESICS, INC. KNOWLEDGE CENTER'. Below this is the heading 'Your Resource for Training and Technical Support'. A search results section titled 'Search results for:' shows a search bar with the text 'location'. Below the search bar, there is a list of search results, each with a bullet point and a title. The results include: 'How to set up a Timing Test for Net Station 5, E-Prime and EENS', 'GPS 3.0 Solver: Video Tutorial', 'How to Uninstall Net Station 4.5.x', 'Spline Interpolation of the Scalp EEG', '3D Sensor Position files for visualizing data in third party software', 'Relating the HC GSN sensor positions to the 10-10 International electrode placement system', 'How to de-identify a file in Net Station', and 'GeoSource 3.0 Research Basic: Video Tutorial'. Each result is followed by a brief description of the content.

EGI ELECTRICAL GEODESICS, INC. KNOWLEDGE CENTER

Your Resource for Training and Technical Support


Search results for:

location


- [How to set up a Timing Test for Net Station 5, E-Prime and EENS](#)
In the assembled photocell holder. g. Position the photocell so the tip is touching the subject's monitor in the location where the stimulus will appear on the screen. 2. Configuring the DIN interface in NetStation 5.3 Acquisition a. Launch Net Station...
- [GPS 3.0 Solver: Video Tutorial](#)
Once you have acquired a GPSR image file using the GPS 3.0 Acquisition System, solving the image file to derive the individual sensor coordinate information is the next step. This video tutorial will walk you through all the steps you need to know to...
- [How to Uninstall Net Station 4.5.x](#)
after reinstall. To save your EEG files, ensure that you move all files within the Net Station User Data folder to another location on your computer. First, navigate to the following locations and remove the following items: a. Macintosh HD/Applications...
- [Spline Interpolation of the Scalp EEG](#)
electrode. To render the measured potentials graphically, a method for estimating scalp potentials in between electrode locations is needed. As described in EGI's Technical Note "Spline Interpolation of the Scalp EEG", the method was first developed in...
- [3D Sensor Position files for visualizing data in third party software](#)
Geodesic Sensor Net™ (HC GSN) into the destination software. The 3D sensor coordinates are stored in several different locations so that both Net Station™ and GeoSource™ 3.0 Research can make use of them. Within Net Station, they are stored in .xml...
- [Relating the HC GSN sensor positions to the 10-10 International electrode placement system](#)
The 10-20 system is an internationally recognized method that describes EEG electrode locations on the head surface in terms of inter-electrode distances relative to distances between cranial landmarks over the head surface. The primary purpose of the...
- [How to de-identify a file in Net Station](#)
to change the name as necessary. Click "Return" when you are finished. *Note that all edits you make across the different locations must be consistent, i.e. if you are changing "Doe,John" to "JD" then this must be the same across all files indicated in...
- [GeoSource 3.0 Research Basic: Video Tutorial](#)

Then choose the file you need to use. Download it and put it on your desktop. Then we need to load it into your EEGLAB data.

In eeglab, go to Edit and click "channel locations", click cancel button in the following pop up window. "Edit channel info window" will appear then.




Then Click “Read Locations” button and choose the location file you just downloaded. Choose the 5th option, so that the .sfp file can be recognized.




Go to the last channel. Then you will see it should be Cz. Uncheck the option “Channel in data array”, and click the “Set reference” button. In the following window: you need to enter “1:257” into Channel indices, and “Cz” into Reference.

Note, Here we use a 256ch data, if your data use different channels, it need to be different. For example, if your data is 128ch, enter “1:129”


After that, you will see your data has be recognized to have ref channel: Cz




Step 5: Average reference(optional). Tools/Re-reference:



The important thing you need to keep in mind is that you need to add ref channel back to your data channel. Check “compute average reference” and add Cz back.



Then you will see that you already change your reference to Average reference. If you want to use another reference, just repeat this step and choose another reference channel whatever you want.



Topic 3 Bad channel replacement

Sometimes we want to make use of ICA function to remove artifacts such as eyeblink, line noise, EKG etc, while we found some bad channels and want to make use of the “bad channel replacement tool” in Netstation to maintain flexibility. Then we do the following steps:

Step 1: filter your data.

This step is necessary since we will use “artifact detection” tool to detect the bad channels in the segments. If you do not remove your DC signals, then this tool cannot work well.

Note: High pass needed. Low pass is optional, you can make use of the filter in EEGLAB or do ICA analysis to remove the line noise.

Net Station Tools

Idle

Tool Settings Jobs/Results Log

Filtering Specification

Specification Name: Highpass

Save

Output Options

Name: Append Operation Name: _fil

Destination: Same As Source

☒ Include video if appropriate

Duplicate

Add to Set

Filtering Settings:

☐ Auto set name: 0.10 Hz Highpass

☒ Highpass

☐ Lowpass


☐ Notch

Options


0.10

Tool Settings


Step 2: Segment. We do segmentation here:



Step 3: Mark up events. As we described in the topic 1, we change the names of the stimulus into different names so that we can analyze them better in EEGLAB




Step 4-artifact detection. We only do bad channel detection here.



Step5 bad channel replacement. Since you do not need to set any parameters to this tool, just create one and run it to the file produced by Step 4.

Step 6. Then repeat the steps in Topic 2. The only difference here is the settings of the option of "file export" tools. We do not need to export the bad segments which are recognized by Nestation: Check the option "export only good segments"



Topic 4 ICA analysis in EEGLAB

Sometimes the customer want to make use of the ICA function in EEGLAB, for example, to remove the artifacts. Netstaion will give up a epoch when there is an Eye blink or Eye movement in it. However, in some experiments, the trail will be last very long(several seconds even longer). So they just want to retain the epoch.

After we follow the step in Topic 2 or Topic 3, we can do ICA analysis in EEGLAB. Here we do not teach how to use ICA in EEGLAB, please refer to the EEGLAB manual which can be downloaded from the official website. We introduce how to run binica to EGI datas.

There are 4 ICA methods in EEGLAB. The default one is runica. Binica is runica's binary code, but can run much faster than Runica.

Step 1 download Binica from the website.

Go to <https://sccn.ucsd.edu/wiki/Binica> .

From there, you can download the most recent packages. Since EGI use Mac as its workstation, we only introduce how to run Binica on Mac. According to CPU type, download the correct one:

Page [Discussion](#) [Read](#) [View source](#) [View history](#)

Binica

Download BINICA binaries: a Compiled C Version of runica()

This binary version of the `runica()` function of Makeig et al. contained in the EEG/ICA Toolbox runs ~12x faster than the Matlab version and may be ~4x more compact. It uses the logistic infomax ICA algorithm of Bell and Sejnowski, with natural gradient and extended ICA extensions. It was programmed for unsupervised operation by Scott Makeig at CNL, Salk Institute, La Jolla CA. Sigurd Enghoff translated it into C++ code and compiled it for multiple platforms. J-R Duann has improved the PCA dimension-reduction and has compiled the linux and free_bsd versions.

Use: To use the function, call with a ".sc" file argument. For individual applications, copy and modify the sample script "ica.sc".m. Under UNIX or DOS command line

```
% ica < myversion.sc
```

Be sure that the directory that you store the binary file in is in your search path. e.g., In Unix/Linux add this directory to your root .cshrc file "setenv path" line.

Outputs: Ica creates two files, "xxx.wts" and "xxx.sph" containing weights and sphere matrices such that (under Matlab)

```
>> ICA_activations = wts * sph * data;
```

The "xxx" stem in the output files may be specified within the input .sc parameter file. See the sample .sc file for arguments, and the EEG/ICA toolbox tutorial for more details.

- [Linux version \(works under Red Hat 7 or 8\)](#) ☐. See also OSX compiled version with GCC below.
- [FreeBSD 4.0](#) ☐
- [Mac OSX PPC](#) ☐ (Recompiled by [William Beaudot](#) ☐ in 2004) or [Recent compile \(April 2016\)](#)
- [Mac OSX Intel](#) (Recompiled by [Grega Repovš](#) ☐ in 2009)
- Windows: the source code should compile on windows although we do not have a binary available. If you have one, please send it to us.

To use one of these programs from within Matlab (and EEGLAB)

1. download the file and place them in the eeglab directory (you may create a subfolder for them or uncompress them in the function subfolder).
2. edit the `icadefs.m` ☐ file to specify the file name of the executable you intend to use (variable ICABINARY).
3. add the path to the binary function both to your Matlab path and to your Unix path, otherwise the system will not be able to locate the executable file.
4. From the command line, you may use the `binica.m` ☐ function that will call the binary executable. From the EEGLAB graphical interface, run ICA using the 'binica' option of the Tools > Run ICA graphic interface (see the tutorial for how to [compute ICA components](#)).

Older version of the ICA binary are available below (these versions are not compatible with the Matlab `binica.m` ☐ function and **cannot be used directly in Matlab or from EEGLAB** (although see [bug 1604](#) ☐); The ica.sc text configuration file ([sample here](#) ☐) in the archive must be edited manually). [SGI Unix \(older version\)](#) ☐ [Sun Unix \(older version\)](#) ☐ [Windows PC \(95, 98, NT, 2000, XP?, older version\)](#) ☐

Download SOURCE code


The [binica.zip](#) ☐ (13Mb) contains the source code, and [binica_full.zip](#) ☐ (~180Mb) source code plus many binaries. This code is distributed under the GNU GPL license and may not be used for commercial applications. It is copyrighted by the Salk Institute for biological studies and the University of San Diego California. This code can usually compile under most Unix machines.

The binary above for Mac OSX also contains a make file for Mac OSX. Some recommendation below:

1. if you uncompress using winzip, deactivate the "tar smart CR/LF" option in winzip in the menu Option > configuration tab Miscellaneous
2. recompile BLAS (folder CLABPACK\BLAS)
3. recompile LABPACK (CLABPACK folder)
For 2 and 3 it is actually better if you find on the Internet the latest versions of these libraries
4. make the ICA binary file by using the makefile in the main directory
5. Modify the `icadefs.m` Matlab file under EEGLAB so that it points to your binary (in case you want to call it from Matlab).
6. For credits, please quote "binary Infomax ICA by Sigurd Enghoff, based on the Matlab version of Scott Makeig and collaborators. Makeig S, Anthony J. Bell, Tzyy-Ping Jung and Terrence J. Sejnowski, Independent component analysis of electroencephalographic data In: D. Touretzky, M. Mozer and M. Hasselmo (Eds). Advances in Neural Information Processing Svstems 8:145-151 (1996)."

Step 2 install Binica to your EEGLAB

Uncompress the zip file to the functions folder under EEGLAB folders. Then copy the “ica_osx” file to resource folder.



Step 3 : go to sigprofunc folder and fine icadefs.m,modify the following sentence

from

```
ICABINARY = fullfile(eeglab_p, 'functions', 'resources', 'ica_linux');
```

To

```
ICABINARY = fullfile(eeglab_p, 'functions', 'resources', 'ica_osx');
```

Save and close.